

Anti-RNA Polymerase III (RDL)

Test Details

METHODOLOGY

Enzyme immunoassay (EIA)

RESULT TURNAROUND TIME

14 - 21 days

Turnaround time is defined as the usual number of days from the date of pickup of a specimen for testing to when the result is released to the ordering provider. In some cases, additional time should be allowed for additional confirmatory or additional reflex tests. Testing schedules may vary.

Use

The detection of anti-RNAP III antibodies is useful in the diagnosis of SSc and for the identification of patients at risk for developing progressive skin thickening and renal crisis. The prevalence of IgG RNAP III antibodies is from 3 - 58% in SSc patients.

Specimen Requirements

Specimen

Serum

Volume

1 mL

Minimum Volume

0.3 mL (**Note:** This volume does **not** allow for repeat testing.)

Container

Red-top tube or gel-barrier tube

Collection Instructions

Separate serum from cells within one hour of collection. Transfer to a plastic transport tube before shipping. To avoid delays in turnaround time when requesting multiple tests on frozen samples, **please submit separate frozen specimens for each test requested.**

Stability Requirements

Temperature	Period
Room temperature	7 days
Refrigerated	14 days
Frozen	60 days
Freeze/thaw cycles	Stable x1

Reference Range

- Negative: <20 Units
 - Weak Positive: 20–39 Units
 - Moderate Positive: 40–80 Units
 - Strong Positive: >80 Units
-

Storage Instructions

Refrigerate or **freeze**.

Causes for Rejection

Grossly hemolyzed; bacterial contamination; lipemic specimen; icteric specimen; non-serum specimen types
