

Apolipoprotein A-1

Test Details

METHODOLOGY

Immunologic

RESULT TURNAROUND TIME

1 - 2 days

Turnaround time is defined as the usual number of days from the date of pickup of a specimen for testing to when the result is released to the ordering provider. In some cases, additional time should be allowed for additional confirmatory or additional reflex tests. Testing schedules may vary.

Use

Apo A-1 comprises 60% to 70% of the protein content of HDL. Apo A-1 serves as an activator of LCAT which is involved in the transport of cholesterol from peripheral tissues to the liver where it is degraded. Decreased serum HDL cholesterol levels have been reported to correlate with increased risk of coronary artery disease (CAD); however, Apo A-1 has been suggested as a better discrimination of CAD than HDL.¹ Apo A-1 levels also correlate with survival rates or risk factors for patients with myocardial infarction and peripheral vascular disease. Measurement of Apo A-1 may also aid in the diagnosis of Tangier disease (absence of alpha-lipoprotein).

Special Instructions

State patient's sex on the test request form.

Specimen Requirements

Specimen

Serum (preferred) **or** plasma

Volume

2 mL

Container

Red-top tube, gel-barrier tube, green-top (heparin) tube, **or** lavender-top (EDTA) tube

Collection Instructions

Separate serum or plasma from cells. Transfer specimen to a plastic transport tube.

Stability Requirements

Temperature	Period
Room temperature	14 days
Refrigerated	14 days
Frozen	14 days
Freeze/thaw cycles	Stable x3

Reference Range

See table.

Male (mg/dL)	
0 to 14 d	62–91
15 d to 11 m	53–175
1 to 5 y	80–164
6 to 12 y	107–175
>12 y	101–178
Female (mg/dL)	
0 to 14 d	71–97
15 d to 11 m	53–175
1 to 5 y	80–164
6 to 80 y	116–209
>80 y	114–214

Storage Instructions

Maintain specimen at room temperature.

Patient Preparation

Patient must be fasting 12 to 14 hours.

Causes for Rejection

Specimen from nonfasting patient

References

Colantonio DA, Kyriakopoulou L, Chan MK, et al. Closing the gaps in pediatric laboratory reference intervals: A CALIPER database of 40 biochemical markers in a healthy and multiethnic population of children.

Clin Chem. 2012 May; 58(5):854-868. PubMed 22371482

(<http://www.ncbi.nlm.nih.gov/pubmed/22371482>)

Footnotes

1. Maciejko JJ, Holmes DR, Kottke BA, et al. Apolipoprotein A-1 as a marker of angiographically assessed coronary-artery disease. *N Engl J Med.* 1983; 309(7):385-389. PubMed 6410239

(<http://www.ncbi.nlm.nih.gov/pubmed/6410239>)