

Creatine Kinase (CK) and Lactate Dehydrogenase (LD), Totals and Isoenzymes

Test Details

METHODOLOGY

Total: kinetic; Isoenzymes: agarose gel electrophoresis with densitometry

RESULT TURNAROUND TIME

4 - 7 days

Turnaround time is defined as the usual number of days from the date of pickup of a specimen for testing to when the result is released to the ordering provider. In some cases, additional time should be allowed for additional confirmatory or additional reflex tests. Testing schedules may vary.

Use

Evaluation of possible myocardial damage or other cardiac diseases. See individual test descriptions for CK Isoenzymes and LD Isoenzymes.

Special Instructions

State sex of patient on the test request form.

Specimen Requirements

Specimen

Serum

Volume

3 mL refrigerated serum **and** 3 mL **room temperature** serum

Container

Red-top tube or gel-barrier tube

Storage Instructions

Creatine kinase: refrigerate. Lactate dehydrogenase: maintain specimen at room temperature. Do **not** freeze or refrigerate. Creatine kinase: **stable** refrigerated for 48 hours or frozen for two weeks. Lactate dehydrogenase: **stable** at room temperature for seven days or refrigerated for no more than three days.

Causes for Rejection

Creatine kinase: hemolysis. Lactate dehydrogenase: hemolysis; prolonged contact of serum with red cells; frozen specimen; specimen refrigerated more than three days after collection at time of testing.
