

Lyme Disease, Line Blot

Test Details

METHODOLOGY

Line blot

RESULT TURNAROUND TIME

1 - 5 days

Turnaround time is defined as the usual number of days from the date of pickup of a specimen for testing to when the result is released to the ordering provider. In some cases, additional time should be allowed for additional confirmatory or additional reflex tests. Testing schedules may vary.

Test Includes

Line blot analysis and interpretations for IgG- and IgM-specific antibodies

Use

Detect antibodies specific for *B burgdorferi*. This test is intended as the second step in testing serum samples found to have been equivocal or positive using an ELISA or IFA test as the first step to provide supportive evidence of infection with *B burgdorferi*.

Special Instructions

This test is intended as the second step in testing serum samples found to have been equivocal or positive using an ELISA or IFA test as the first step to provide supportive evidence of infection with *B burgdorferi*. For CSF, order test **160457** (/tests/160457/lyme-disease-line-blot-cerebrospinal-fluid).

Custom Additional Information

Provides both IgM and IgG Line blot results. Increased specificity reduces cross-reactivity from other spirochetes. Detect and analyze antibodies against specific *B burgdorferi* proteins. The sensitivity of a Line blot assay by itself is not sufficient to assess Lyme infection adequately. Both EIA and Line blot tests should be performed.

Specimen Requirements

Specimen

Serum

Volume

1 mL

Minimum Volume

0.5 mL

Container

Red-top tube or gel-barrier tube or serum transfer tube

Stability Requirements

Temperature	Period
Room temperature	14 days
Refrigerated	14 days
Frozen	14 days
Freeze/thaw cycles	Stable x3

Reference Range

Negative

Storage Instructions

Maintain specimen at room temperature.

Causes for Rejection

Hemolysis; lipemia; gross bacterial contamination
