

Ovarian Cancer Monitor

Test Details

METHODOLOGY

HE4: enzyme-linked immunosorbent assay (ELISA); CA 125: electrochemiluminescence immunoassay (ECLIA)

RESULT TURNAROUND TIME

1 - 3 days

Turnaround time is defined as the usual number of days from the date of pickup of a specimen for testing to when the result is released to the ordering provider. In some cases, additional time should be allowed for additional confirmatory or additional reflex tests. Testing schedules may vary.

Test Includes

Human Epididymis Protein 4 [081700] (/tests/081700/human-epididymis-protein-4); Cancer Antigen 125 [002303] (/tests/002303/cancer-antigen-ca-125). See individual test descriptions for complete information.

Use

Use of HE4 in combination with CA 125 provides a more accurate predictor of malignancy than either marker alone.

Special Instructions

This test may exhibit interference when sample is collected from a person who is consuming a supplement with a high dose of biotin (also termed as vitamin B7 or B8, vitamin H, or coenzyme R). It is recommended to ask all patients who may be indicated for this test about biotin supplementation. Patients should be cautioned to stop biotin consumption at least 72 hours prior to the collection of a sample.

Specimen Requirements

Specimen

Serum

Volume

1 mL

Minimum Volume

0.8 mL (**Note:** This volume does **not** allow for repeat testing.)

Container

Red-top tube or gel-barrier tube

Collection Instructions

If red-top tube is used, transfer separated serum to a plastic transport tube.

Stability Requirements

Temperature	Period
Room temperature	5 days
Refrigerated	14 days
Frozen	14 days
Freeze/thaw cycles	Stable x3

Storage Instructions

Refrigerate
