

Syphilis: RPR With Reflex to RPR Titer

Test Details

METHODOLOGY

Charcoal flocculation

RESULT TURNAROUND TIME

Within 1 day

Turnaround time is defined as the usual number of days from the date of pickup of a specimen for testing to when the result is released to the ordering provider. In some cases, additional time should be allowed for additional confirmatory or additional reflex tests. Testing schedules may vary.

Use

This test can be used to follow treatment response in patients being treated for syphilis infection. This assay should not be used to screen for syphilis infection. To screen for syphilis infection, a reflex cascade that includes both RPR and a *treponema*-specific assay should be utilized, such as *Treponema pallidum* (Syphilis) Screening Cascade [082345] or Rapid Plasma Reagin (RPR) Test With Reflex to Quantitative RPR and Confirmatory *Treponema pallidum* Antibodies [012005]. See Related Information to view these tests.

Specimen Requirements

Specimen

Serum

Volume

1 mL

Minimum Volume

0.5 mL

Container

Red-top tube or gel-barrier tube

Stability Requirements

Temperature	Period
Room temperature	14 days
Refrigerated	14 days
Frozen	14 days
Freeze/thaw cycles	Stable x3

Reference Range

Nonreactive

Storage Instructions

Room temperature

Causes for Rejection

Hemolysis; lipemia; gross bacterial contamination; plasma specimen
