

Paramount Diagnostic Lab

Test Details

METHODOLOGY

Electrochemiluminescence immunoassay (ECLIA)

RESULT TURNAROUND TIME

Within 1 day

Turnaround time is defined as the usual number of days from the date of pickup of a specimen for testing to when the result is released to the ordering provider. In some cases, additional time should be allowed for additional confirmatory or additional reflex tests. Testing schedules may vary.

Test Includes

Folate; vitamin B₁₂

Special Instructions

This test may exhibit interference when sample is collected from a person who is consuming a supplement with a high dose of biotin (also termed as vitamin B7 or B8, vitamin H, or coenzyme R). It is recommended to ask all patients who may be indicated for this test about biotin supplementation. Patients should be cautioned to stop biotin consumption at least 72 hours prior to the collection of a sample.

Specimen Requirements

Specimen

Serum

Volume

1.5 mL

Minimum Volume

0.6 mL (**Note:** This volume does **not** allow for repeat testing.)

Container

Red-top tube or gel-barrier tube

Collection Instructions

If red-top tube is used, transfer separated serum **immediately** to a plastic transport tube.

Stability Requirements

Temperature	Period
Room temperature	72 hours
Refrigerated	7 days
Frozen	14 days
Freeze/thaw cycles	Stable x3

Storage Instructions

Room temperature or refrigerate

Causes for Rejection

Plasma specimen
