

ANA 12 Plus Profile (RDL)

Test Details

METHODOLOGY

ELISA, IFA, Farr RIA, Turbidimetry, CMIA

RESULT TURNAROUND TIME

7 - 10 days

Turnaround time is defined as the usual number of days from the date of pickup of a specimen for testing to when the result is released to the ordering provider. In some cases, additional time should be allowed for additional confirmatory or additional reflex tests. Testing schedules may vary.

Test Includes

ANA 12 Plus Profile is reflexive and consists of the following: Anti-Nuclear Ab (ANA). If positive, reflexes to following 14 components: Anti-Cardiolipin Ab (IgG, IgA & IgM Isotypes); Anti-Centromere Ab; Anti-Chromatin Ab; Anti-dsDNA Ab by Farr; Anti-La (SS-B) Ab; Anti-Ro (SS-A) Ab; Anti-Scl-70 Ab; Anti-Sm; Anti-TPO (Thyroid Microsomal Peroxidase) Ab; Anti-U1 RNP; C3 Complement; C4 Complement; Rheumatoid Factor by Turb; Anti-CCP (Cyclic Citrullinated Peptide) Ab.

If ANA is positive, all tests listed are performed. If ANA is negative, an Anti-Ro (SS-A) Ab test is performed to catch the small percentage of ANA false negative results. ANA includes titer & pattern.

ANA 12 Plus Profile, Do All is not reflexive and is a separate test code (use 520175).

Use

The ANA 12 Plus Profile is a reflexive assay. If the antinuclear antibodies titer is negative, no further tests are performed except Anti-Ro Ab. Anti-Ro is performed to find the very small percentage of false negative antinuclear antibodies. If the antinuclear antibodies titer is positive, the complete profile is performed. The antinuclear antibodies test detects autoantibodies that bind to the nuclear components of HEp-2 cells. A positive antinuclear antibodies titer is not specific but indicates the presence of one or more autoantibodies. Profiles of specific autoantibodies have diagnostic and potential prognostic utility.

This profile helps differentiate between eight autoimmune disease states:

SLE

- Antinuclear antibodies by IFA

Tests with high specificity

- Anti-dsDNA by Farr
- Anti-SM

Tests with moderate specificity

- Anti-Ro
- Anti-La
- Anti-Chromatin
- C3 & C4 Complements
- Anticardiolipin IgG, IgA & IgM

Drug-induced SLE

- Anti-Chromatin

Sjogren's Syndrome

- Anti-Ro
- Anti-La

Rheumatoid Arthritis

- Anti-CCP (high specificity)
- Rheumatoid Factor (Turbidimetry)

Mixed Connective Tissue Disease (MCTD)

- Anti-U1 RNP

Limited Scleroderma (CREST)

- Anti-Centomere

Diffuse Scleroderma

- Anti-SCL-70 by EIA, confirmed by a second method

Autoimmune Thyroid Disease

- Anti-Thyroid Microsomal (TPO)

Special Instructions

If the antinuclear antibodies titer is positive, all tests listed are performed.

If the antinuclear antibodies titer is negative, an Anti-Ro (SS-A) Ab test is performed to catch the small percentage of false negative antinuclear antibodies results.

Antinuclear antibodies includes titer and pattern. ANA 12 Plus is reflexive and consists of the following: Anti-Nuclear Ab (ANA), AntidsDNA Ab by Farr, Anti-ENA Abs (anti-Sm & Anti-U1-RNP), Anti-Centromere Ab, C3 & C4 Complements, Anti-Ro (SS-A) Ab, Anti-La (SS-B), Anti-SCL-70 Ab (IgG, IgA, & IgA isotypes), Anti-TPO (Thyroid Microsomal Peroxidase) Ab, Anti-Chromatin Ab, Rheumatoid Factor by Turb, Anti-CCP (Cyclic Citrullinated Peptide) Ab. To order ANA 12 Plus Profile, Do All, which is not reflexive, please use Test No. 520175.

Specimen Requirements

Specimen

Serum

Volume

5 mL

Minimum Volume

4 mL (**Note:** This volume does **not** allow for repeat testing.)

Container

Red-top tube or gel-barrier tube

Collection Instructions

Separate serum from cells within one hour of collection. Transfer to a plastic transport tube before shipping.

Stability Requirements

Temperature	Period
Room temperature	1 day
Refrigerated	8 days
Frozen	8 days

Storage Instructions

Refrigerate.

Causes for Rejection

Grossly hemolyzed, bacterial contamination, lipemic, icteric, non-serum specimen types
