

Anti-TPO Ab (RDL)

Test Details

METHODOLOGY

Chemiluminescent Immunoassay

RESULT TURNAROUND TIME

14 - 21 days

Turnaround time is defined as the usual number of days from the date of pickup of a specimen for testing to when the result is released to the ordering provider. In some cases, additional time should be allowed for additional confirmatory or additional reflex tests. Testing schedules may vary.

Test Includes

This test is a component of ANA Profile 12, ANA 12 Plus Profile, Thyroid Autoantibody Panel.

Use

In virtually all cases of Hashimoto's disease and in the majority of Graves' disease cases, Anti-TPO autoantibodies are elevated. High levels of Anti-TPO autoantibodies, in the context of the clinical presentation of hypothyroidism, confirms the diagnosis of Hashimoto's disease.

Specimen Requirements

Specimen

Serum

Volume

1 mL

Minimum Volume

0.3 mL (**Note:** This volume does **not** allow for repeat testing.)

Container

Red-top tube or gel-barrier tube

Collection Instructions

Separate serum from cells within one hour of collection. Transfer to a plastic transport tube before shipping. To avoid delays in turnaround time when requesting multiple tests on frozen samples, **please submit separate frozen specimens for each test requested.**

Stability Requirements

Temperature	Period
Room temperature	7 days
Refrigerated	14 days
Frozen	60 days
Freeze/thaw cycles	Stable x1

Storage Instructions

Refrigerate or **freeze**.

Causes for Rejection

Grossly hemolyzed; bacterial contamination; lipemic specimen; icteric specimen; non-serum specimen types
