

Complete Blood Count (CBC) Without Differential

Test Details

METHODOLOGY

Automated cell counter

RESULT TURNAROUND TIME

Within 1 day

Turnaround time is defined as the usual number of days from the date of pickup of a specimen for testing to when the result is released to the ordering provider. In some cases, additional time should be allowed for additional confirmatory or additional reflex tests. Testing schedules may vary.

Test Includes

Hematocrit; hemoglobin; MCH; MCHC; MCV; RDW; red blood cell count (RBC); platelet count; white blood cell count (WBC)

Use

See comments in Complete Blood Count (CBC) With Differential [005009].

Specimen Requirements

Specimen

Whole blood

Volume

Fill tube to capacity.

Minimum Volume

0.5 mL (500 µL for pediatric microtainer capillary tubes; fill tube to capacity)

Container

Lavender-top (EDTA) tube

Collection Instructions

Invert tube 8 to 10 times immediately after tube is filled at the time of collection.

Stability Requirements

Temperature	Period
Room temperature	1 day
Refrigerated	3 days
Frozen	Unstable
Freeze/thaw cycles	Unstable

Storage Instructions

Maintain specimen at room temperature.

Causes for Rejection

Hemolysis; tube not filled with minimum fill volume; specimen drawn in any anticoagulant other than EDTA; specimens diluted or contaminated with IV fluid; clotted specimen; improper labeling; transfer tubes with whole blood; lavender-top (EDTA) tubes received with plasma removed
