

Corticosteroid-binding Globulin (CBG)

Test Details

METHODOLOGY

Radioimmunoassay (RIA)

RESULT TURNAROUND TIME

6 - 10 days

Turnaround time is defined as the usual number of days from the date of pickup of a specimen for testing to when the result is released to the ordering provider. In some cases, additional time should be allowed for additional confirmatory or additional reflex tests. Testing schedules may vary.

Specimen Requirements

Specimen

Serum, **frozen**

Volume

1 mL

Minimum Volume

0.5 mL (**Note:** This volume does **not** allow for repeat testing.)

Container

Red-top tube or gel-barrier tube

Collection Instructions

Separate serum from cells within one hour of collection. Transfer specimen to a plastic transport tube before freezing. To avoid delays in turnaround time when requesting multiple tests on frozen samples, **please submit separate frozen specimens for each test requested.**

Stability Requirements

Temperature	Period
Room temperature	2 days
Refrigerated	2 days
Frozen	300 days

Storage Instructions

Freeze.
