

Human Immunodeficiency Virus 1 (HIV-1) DNA Sequencing Integrase

Test Details

METHODOLOGY

Polymerase chain reaction (PCR) amplification and next generation DNA sequencing

RESULT TURNAROUND TIME

10 - 14 days

Turnaround time is defined as the usual number of days from the date of pickup of a specimen for testing to when the result is released to the ordering provider. In some cases, additional time should be allowed for additional confirmatory or additional reflex tests. Testing schedules may vary.

Use

Detect resistance of HIV-1 to integrase inhibitors in cell-associated viral DNA.

Limitations

This procedure should be used for patients with documented HIV-1 infection and undetectable viral load or low level viremia.

Specimen Requirements

Specimen

Whole blood, **frozen**

Volume

4 mL

Collection Instructions

Collect specimen in lavender-top (EDTA) tube. Do **not** centrifuge. Ship frozen. To avoid delays in turnaround time when requesting multiple tests on frozen samples, **please submit separate frozen specimens for each test requested.**

Storage Instructions

Freeze.

Causes for Rejection

Quantity not sufficient for analysis; heparinized plasma; nonfrozen specimen

References

Günthard HF, Aberg JA, Eron JJ, et al. Antiretroviral treatment of adult HIV infection: 2014 recommendations of the International Antiviral Society-USA Panel. *JAMA*. 2014 Jul 23-30; 312(4):410-425. PubMed 25038359

(<http://www.ncbi.nlm.nih.gov/pubmed/25038359>)

Hurt C. Transmitted resistance to HIV integrase strand-transfer inhibitors: Right on schedule. *Antivir Ther*. 2011; 16(2):137-140. PubMed 21447861

(<http://www.ncbi.nlm.nih.gov/pubmed/21447861>)

Panel on Antiretroviral Guidelines for Adults and Adolescents—A Working Group of the Office of AIDS Research Advisory Council (OARAC). *Guidelines for the Use of Antiretroviral Agents in HIV-1-Infected Adults and Adolescents*.

(<http://www.aidsinfo.nih.gov/ContentFiles/AdultandAdolescentGL.pdf>) Washington, DC:
Department of Health and Human Services; May 1, 2014.