

Paramount Diagnostic Lab

Iron

Test Details

METHODOLOGY

Colorimetric

RESULT TURNAROUND TIME

Within 1 day

Turnaround time is defined as the usual number of days from the date of pickup of a specimen for testing to when the result is released to the ordering provider. In some cases, additional time should be allowed for additional confirmatory or additional reflex tests. Testing schedules may vary.

Use

Aid in the evaluation of a number of conditions involving red cell production and destruction, iron metabolism, or iron transport

Specimen Requirements

Specimen

Serum (preferred) **or** plasma (lithium heparin only). EDTA plasma causes decreased results.

Volume

1 mL

Minimum Volume

0.7 mL (**Note:** This volume does **not** allow for repeat testing.)

Container

Red-top tube, gel-barrier tube, **or** green-top (lithium heparin) tube. Do **not** use oxalate, EDTA, or citrate plasma.

Collection Instructions

Separate serum or plasma from cells within 45 minutes of collection.

Stability Requirements

Temperature	Period
Room temperature	14 days
Refrigerated	14 days
Frozen	14 days
Freeze/thaw cycles	Stable x3

Reference Range

See table.

Male	Range (µg/dL)
0 to 30 d	35–160
1 m to 1 y	18–126
2 to 12 y	28–147
13 to 17 y	26–169
>17 y	38–169
Female	Range (µg/dL)
0 to 30 d	27–133
1 m to 1 y	18–126
2 to 12 y	28–147
13 to 17 y	26–169
18 to 60 y	27–159
>60 y	27–139

Storage Instructions

Maintain specimen at room temperature.

Causes for Rejection

Severe hemolysis; improper labeling; improper collection; EDTA or oxalate plasma specimens
