

Respiratory Pathogen Panel

Test Details

METHODOLOGY

Polymerase chain reaction (PCR)

RESULT TURNAROUND TIME

2 - 4 days

Turnaround time is defined as the usual number of days from the date of pickup of a specimen for testing to when the result is released to the ordering provider. In some cases, additional time should be allowed for additional confirmatory or additional reflex tests. Testing schedules may vary.

Test Includes

Adenovirus; *Bordetella pertussis*; *Bordetella parapertussis*; *Chlamydophila pneumoniae*; SARS-CoV-2; coronavirus 229E; coronavirus HKU1; coronavirus NL63; coronavirus OC43; human metapneumovirus; human rhinovirus/enterovirus; influenza A; influenza A subtype H1; influenza A subtype H3; influenza A subtype H1-2009; influenza B; *Mycoplasma pneumoniae*; parainfluenza virus 1; parainfluenza virus 2; parainfluenza virus 3; parainfluenza virus 4; respiratory syncytial virus

Use

This test provides simultaneous qualitative detection and identification of multiple respiratory viral and bacterial nucleic acids in nasopharyngeal swabs obtained from individuals suspected of respiratory tract infections.

Note: Labcorp also offers smaller targeted panels for Influenza, RSV and SARS-CoV-2 (see test numbers 140140, 140147 or 140163).

Limitations

The performance of this test has not been established for patients without signs and symptoms of respiratory infection. Results from this test must be correlated with the clinical history, epidemiological data, and other data available to the clinician evaluating the patient. Viral and bacterial nucleic acids may persist in vivo independent of organism viability.

Detection of organism target(s) does not imply that the corresponding organisms are infectious or are the causative agents for clinical symptoms. The detection of viral and bacterial nucleic acid is dependent upon proper specimen collection, handling, transportation, storage and preparation. Failure to observe proper procedures in any one of these steps can lead to incorrect results. There is a risk of false positive or false negative values resulting from improperly collected, transported or handled specimens.

Specimen Requirements

Specimen

Nasopharyngeal swab

Volume

One swab in viral or universal transport (VTM/UTM) medium

Minimum Volume

One swab in 0.3 mL transport medium

Container

Viral transport medium

Collection Instructions

Collect NP swab per standard techniques and immediately place swab in up to 3 mL of viral transport media.

Stability Requirements

Temperature	Period
Room temperature	4 hours (stability provided by manufacturer or literature reference)
Refrigerated	3 days (stability provided by manufacturer or literature reference)
Frozen	30 days at 15°C or below (stability provided by manufacturer or literature reference)

Storage Instructions

Freeze.

Causes for Rejection

Improper collection; inadequate specimen; improper labeling; specimen leaked in transit; quantity not sufficient for analysis; name discrepancies; specimen exceeding storage requirements; nasal swabs or any other upper respiratory sample other than nasopharyngeal swab in viral transport media
