

Varicella Zoster Virus (VZV) Antibodies, IgG

Test Details

METHODOLOGY

Chemiluminescent immunoassay (CLIA)

RESULT TURNAROUND TIME

1 - 2 days

Turnaround time is defined as the usual number of days from the date of pickup of a specimen for testing to when the result is released to the ordering provider. In some cases, additional time should be allowed for additional confirmatory or additional reflex tests. Testing schedules may vary.

Use

This test is used for qualitative detection of specific IgG antibodies to varicella-zoster virus (VZV) to aid in the identification of previous infection with VZV and/or vaccination against VZV.

Specimen Requirements

Specimen

Serum

Volume

0.5 mL

Minimum Volume

0.2 mL

Container

Red-top tube or gel-barrier tube

Collection Instructions

If a tube other than a gel-barrier tube is used, transfer separated serum to a plastic transport tube. Do **not** freeze gel-barrier tube (pour off serum first).

Stability Requirements

Temperature	Period
Room temperature	14 days
Refrigerated	14 days
Frozen	14 days
Freeze/thaw cycles	Stable x4

Reference Range

Nonreactive

Storage Instructions

Room temperature

Causes for Rejection

Hemolysis; lipemia; gross bacterial contamination
